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# Rhinogen<sup>®</sup> Bio-Bright<sup>™</sup> One-Step 萤火虫荧光素酶分析试剂盒

货号：RA-GL07

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---

1. De Wet, J.R. et al. (1987) Firefly luciferase gene: Structure and expression in mammalian cells. *Mol. Cell. Biol.* 7, 725–37.
  2. Tatsumi, H., et al. (1989). Luciferase cDNA from Japanese firefly, *Luciola cruciata*: cloning, structure and expression in *Escherichia coli*. *J Biolumin Chemilumin* 3 (2) :75-8.
  3. Collin Goddard (1994) . Cell based screening approaches: advantages of highly automated robotics technology in HTS. *Handbook for the 1994 International Forum on Advances in Screening Technologies and Data Management*, p.19.
  4. Wood, K.V. (1990) Firefly luciferase: A new tool for molecular biologists. *Promega Notes* 28, 1–3.
  5. Alam, J. and Cook, J.L. (1990) Reporter genes: Application to the study of mammalian gene transcription. *Anal. Biochem.* 188, 245–54.
  6. Fraga, H. (2008) . Firefly luminescence: A historical perspective and recent developments. *Photochemical & Photobiological Science*, 7 (2) , 146-158.
  7. Parekh, B.S. et al. (2012) Development and validation of an antibody-dependent cell-mediated cytotoxicity reporter gene assay. *mAbs* 4, 310–8.
  8. Cheng, Z.J. et al. (2012) Development of a bioluminescent cell-based bioassay to measure Fc receptor functionality in antibody-dependent cell-mediated cytotoxicity. American Association of Cancer Research (AACR) Annual Meeting, poster # 2840.
  9. Wood., K.V. (1991) In: *Bioluminescence and Chemiluminescence: Current Status*, Stanley, P., and Kricka, L., eds., John Wiley and Sons, Chichester, NY 543.
  10. Ow, D.W. et al. (1986) Transient and stable expression of the firefly luciferase gene in plant cells and transgenic plants. *Science* 234, 856–9.
  11. Bakhtiarova, A. et al. (2006) Resveratrol inhibits firefly luciferase. *Biochem. Biophys. Res. Comm.* 351, 481–4.
-

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