

Certificate of Analysis

Recombinant human Follicle Stimulating Hormone

Cat. No.: QRE-201
Lot. No.: AG2010220301
有效期至: 10 Mar 2025

Storage: The product is shipped with blue ice and should be stored below -20°C upon arrival.

Specificity by ELISA: Conforms.

Protein Purity Assay (SDS-PAGE) : no less than 98% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

Protein Concentration by E280: 1.31 mg/ml

Extinction Coefficient at 0.1% = 1.026

For in vitro biological function research and clinical diagnostic raw materials. Not for Human testing.

Shanghai RhinoBio Co.,Ltd

99 Ruoshui Road, Suzhou Industrial Park, Jiangsu Province .

For inquiries, contact us at 0086-512-87663137or techserv@rhinobio.com.

www.rhinobio.com


QA/QC Manager
Issued on 11 Mar 2022
质检专用章