

# Certificate of Analysis

## Recombinant human Follicle Stimulating Hormone

Cat. No.: QRE-201  
Lot. No.: AG2010230903  
Expiration Date: 30 May 2024

---

**Storage:** The product is shipped with blue ice and should be stored below -70°C upon arrival.

**Specificity by ELISA:** Conforms.

**Protein Purity Assay (SDS-PAGE) :** no less than 98% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Protein Concentration by E280:** 1.18 mg/ml

**Extinction Coefficient** at 0.1% = 1.026

---

For in vitro biological function research and clinical diagnostic raw materials. Not for Human testing.

### Shanghai RhinoBio Co.,Ltd

99 Ruoshui Road, Suzhou Industrial Park, Jiangsu Province .  
For inquiries, contact us at 0086-512-87663137or techserv@rhinobio.com.  
www.rhinobio.com

