## ertificate of Analysis



**Product Name** 

Human FcγR IIIa 158V(His 标签)

Cat. No.

**QRE-137** 

Lot. No.

AG370240501

Storage

The product is shipped with blue ice and should be stored below -20°C upon arrival.

TEST	SPECIFICATIONS	RESULTS LOT
Specificity by ELISA		Pass
Protein Purity Assay (SDS-PAGE)	≥95%	Pass
Protein Concentration by E280	1.0 mg/ml	
Extinction Coefficient at 0.1%	1.849	
Date of manufacture		Apr 2024
Expiration date		Mar 2025



Contact: Phone 0512-87663137

Sales email: sales@rhinobio.com

Official website: www.rhinobio.com

Shanghai Rhinobio Biotechnology Co, Ltd

Building B6, Suzhou Nanobiotechnology Park, Suzhou

Industrial Park